



STAT KIT® Z-1000 Mobile ACLS EMERGENCY MEDICAL KIT

ACLS medications and devices are organized for quick access. Includes pediatric epinephrine auto-injector. Add oxygen and an AED for an integrated, portable emergency response kit.



AED and Oxygen
Sold Separately



STAT KIT® Z-1000 Mobile ACLS highlights

The STAT KIT Z-1000 Mobile ACLS organizes its medications and medical devices into sections, compartments and trays, all clearly labeled and categorized into key areas of care. This allows the responder to immediately assess and treat the patient without searching for devices and medications.

A complete set of medications, devices, oxygen system*, and defibrillation* to provide:

- Basic Life Support
- Advanced Cardiac Life Support (ACLS) for up to 30 minutes



Always up to date

HealthFirst's Auto-Replenishment Service takes the guesswork out of staying up to date. The service uses smart automation to track the expiration dates of medications and devices and automatically replenishes them before expiration. Your practice benefits by extending the shelf life of medications, lowering costs, saving time, and gaining peace of mind.

OnTraq

OnTraq is an online account management tool that dynamically shows the status of your emergency medical kit and medications. It allows users to view each medication, expiration dates, lot numbers, national drug backorder status, and upcoming medication replenishment shipments. Rely on OnTraq's smart automation for dashboards, reports and task management.

The STAT KIT Z-1000 Mobile ACLS contains the following essential medications and devices:

Medications

Anaphylaxis, Allergy & Asthma Tray

Albuterol inhaler 90mcg
Diphenhydramine 25mg capsule (2)
Diphenhydramine 50mg/mL, 2mL
Epinephrine auto-injector (Pedi) 0.15mg
Epinephrine auto-injector (Adult) 0.3mg
Epinephrine 1:1,000 1mg/mL, 2mL
Solu-medrol 125mg, 2mL

Cardiac Tray

Adenosine 3mg/mL, 10mL
Amiodarone 50mg/mL, 6mL
Aspirin tablets 325mg (4)
Atropine sulfate 0.1mg/mL, 20mL
Dopamine 40mg/mL, 10mL
Epinephrine 1:10,000, 50mL
Lidocaine 2% 20mg/mL, 10mL
Magnesium sulfate USP 50% 5g/10mL, 10mL
Nitroglycerin 0.4mg tablets
Verapamil 2.5mg/mL, 4mL

Miscellaneous Drug Tray

Dextrose 25% (Pedi) 250mg/mL, 10mL
Dextrose 50% (Adult) 500mg/mL, 50mL
Flumazenil 0.1mg/mL, 10mL
Naloxone 0.4mg/mL, 2mL
Ondansetron 2mg/mL, 4mL
Oral glucose gel

IV Equipment & Syringes

Syringes

3cc 21G x 1 1/2" (1)
3cc 20G x 1 1/2" (1)
10cc 21G x 1 1/2" (1)
10cc 20G x 1 1/2" (1)
1cc syringe (4)
3cc syringe (1)
12cc syringe (1)
Set of syringe needles

Sodium Chloride 0.9%

500mL bag (2)
IV administration tubing
IV catheter, 20G x 1" (2)
IV catheter, 22G x 1" (1)

Miscellaneous

Alcohol wipes (6)
Nitrile gloves, pair (1)
Tape, hypoallergenic (1)
Tourniquet, latex free (2)
Security seals (4)

Airway Equipment

Oropharyngeal airway set
60mm (1)
80mm (1)
100mm (1)
Laryngeal Mask Airways (LMA)
Size 3 (1)
Size 4 (1)
Size 5 (1)
Syringe for inflation (1)
Endotracheal tube set
3mm (1)
5mm (1)
7mm (1)
Laryngoscope (1)
Size 2 blade (1)
Size 4 blade (1)
Battery, C cell (1 pair)
Hand-held suction unit (1)
Yankauer suction tip (1)
Magill forceps (1)

Resuscitation & Oxygen

Bag-Valve-Mask (BVM) resuscitator (1)
Adult mask (1)
Pediatric mask (1)
CO₂ detector

*Oxygen and AEDs are available.
Ask about adding them to complete your kit.

Trauma

Gauze pads, 4"x 4" (4)
Gauze bandage, rolled (2)
Compression bandage, 3" (2)
Swabstick and antiseptic (1)
Gauze sponges, 3"x 3" (2)
Nitrile gloves, pair (1)
Scalpel (1)
Shears (1)
Cold pack (1)

Vitals, PPE, Guides

Stethoscope (1)
Sphygmomanometer (1)
Thermometer strips (2)
Headlamp (1)
Nitrile gloves, pair (1)
N95 respirator (1)
CPR shield (1)
Safety glasses (1)
ACLS pocket guide
First aid guide
Notepad and pen

A physician license is required for the purchase of prescription medications and devices.

IMPORTANT: Medications will be shipped in the available form factor at the time of fulfillment. Some medications or devices may not be available at the time of order due to national supply issues. Kits may ship with backordered items. Backordered items will ship as soon as they are available. Visit healthfirst.com for complete terms and conditions of sale.