



CONFIRM™ 10 IN-OFFICE SPORE TESTING SYSTEM

Safe. Simple. Compliant.

Final Results in Just 10 Hours

Ensure your sterilization cycles are effective with the Confirm 10 Steam Biological Indicators and reusable incubator system.

- 25 indicators per box
- Detect spore growth in as little as 2.5 hours
- Results visible in 10 hours—yellow indicates a failed test

How It Works:

- 1 Run the indicator through a normal sterilizer cycle
- 2 Activate and incubate
- 3 Record results and re-test if needed

Use OnTraq, HealthFirst's online compliance dashboard, to log and manage your test results with ease.



Also Available from HealthFirst:

In-Office Infection Control and Waste Management Solutions

Track It All in OnTraq™ Your Central Compliance Dashboard

Waterline Testing Solutions

Maintain safe dental unit water with easy-to-use testing kits. Log results and monitor compliance directly in OnTraq.

Sharps Disposal

Protect your team and stay compliant with sharps containers and mail-back disposal services fully trackable in OnTraq.

Pharmaceutical Waste Management

Safely dispose of unused or partially used medications, including dental carpules, with HealthFirst's compliant solutions. Document disposal activity in OnTraq for complete reporting.

Information Technology Built to Protect Patients

compliance management
replenishes medications and devices
on and reporting for easy audits
and readiness
patients

Learn more or order today at healthfirst.com/spore 800-331-1984

11433 v001 01/2026 All orders fulfilled by HF Acquisition Co. LLC