

AED Comparison Chart

	HeartSine		Defibtech	ZOLL	
	Samaritan® PAD 350P/360P	Samaritan® PAD 450P	Lifeline AED®	AED Plus®	AED Pro®
<i>HealthFirst is an authorized distributor of the following AEDs. This chart is for general comparison. Please refer to the manufacturer spec sheets for more details.</i>					
Fully or Semi-Automatic	350P Semi-automatic, 360P Fully automatic	Semi-automatic	Semi-automatic	Semi-automatic	Semi-automatic with manual override
Energy Output (Joules) Adult Pads	150J-150J-200J	150J-150J-200J	150J-150J-150J	120J-150J-200J	120J-150J-200J
Energy Output (Joules) Pediatric Pads	50J	50J	50J	50J, 70J, 85J	50J-70J-85J
Waveform	Biphasic	Biphasic	Biphasic	Biphasic	Biphasic
Pediatric Delivery Method	Pediatric pads (sold separately)	Pediatric pads (sold separately)	Pediatric pads (sold separately)	Pediatric pads (sold separately)	Pediatric pads (sold separately)
AED Prompt/Instruction Type	Voice, LED	Voice, LED	Voice, LED	Voice, LED	Voice, LED, text
CPR Feedback/Coaching	Yes, visual and voice prompts	Yes, visual and voice prompts	Yes, metronome	Yes, visual and voice prompts	Yes, visual and voice prompts
ECG Display	No	No	No	No	Yes, 3.02" x 2.27" LCD display*
ECG Cable	No	No	No	No	3-lead
Electrode Pads Life (years)	4	4	2	5	5
Battery Standby Life (years)	4	4	5, or 7 optional	5	5
Self-test	Yes	Yes	Yes	Yes	Yes
Dimensions	8.0" x 7.25" x 1.9"	8.0" x 7.25" x 1.9"	11.8" x 8.5" x 2.7"	5.25" x 9.50" x 11.50"	9.4" x 9.2" x 3.0"
Weight (pounds)	2.4	2.4	4.2	6.7	6.5
Design/Use	General public	General public	General public	General public	Professional
Data Download Method	USB	USB	Data card	IrDA dongle	Wireless or USB
Warranty (years)	8	8	8	5, 7 with registration	5, 7 with registration
Dust/Water Resistance	IP56	IP56	IP54	IP55	IP55
Initial purchase includes	Adult Pad-Pak with integrated battery	Adult Pak-Pak with integrated battery	Adult pads and 5-year battery	Adult CPR-D PADZ and battery	Pads and battery sold separately
Pads, batteries and AED accessories available:	Batteries, pediatric pads, case, data cable, cabinet, wall bracket, AED trainer, trainer pads, AED wall sign, prep kit	Batteries, pediatric pads, case, data cable, cabinet, wall bracket, AED trainer, trainer pads, AED wall sign, prep kit	Batteries, pediatric pads, data cards, case, AED wall sign, prep kit	Batteries, Zoll CPR-D PADZ, Pediatric pads, mounting brackets, case, wall cabinets, AED wall sign, prep kit	Batteries, Zoll CPR-D PADZ, Pediatric pads, mounting brackets, case, wall cabinets, AED wall sign, prep kit
Manufacturer product downloadable files:	HeartSine PAD 350P/360P	HeartSine PAD 450P	Defibtech Lifeline AED	ZOLL AED Plus	ZOLL AED PRO

*Heart rate, ECG waveform, text prompts, CPR bar graph, battery gauge, elapsed time, number of shocks delivered.

HealthFirst manages your AED pads and batteries, sending replacements automatically before they expire through OnTraq.™

The information contained in this sales sheet is provided for informational purposes only and is subject to change without notice. All products are sold as-is, and the distributor makes no warranties or representations, either express or implied, regarding the accuracy, completeness, or suitability of the products for any particular purpose. All prices and product availability are subject to change. It is the responsibility of the purchaser to verify product specifications and ensure the products are suitable for their intended use.

