

3 Gallon R Style Sharps Container



Product Code: 0330-150R

General Description

- Multi-purpose container suitable for OR and dialysis
- Rotating Door provides Temporary closure
- Locks for final Disposal
- Dual openings accommodates sharps of various sizes.
- Includes needle unwinder
- Puncture resistant
- Made of high impact polypropylene
- Reinforced hinge as requested by FDA
- Meets all FDA, OSHA and NOSHA
- Use instructions on part label
- Comparable to Covidien: 8527R, 8964
- Fits Covidien brackets: 8524C, 8528C, 8529




Rotate door to close. For final disposal rotate door till arrows point at each other.



Packing Details

Part Dimensions:	13.4 x 6 x 13.5
Case:	10 each/case
Packaged:	Bulk, Nested
Case Dimensions:	22.50 x 11.00 x 14.00
Case Volume:	2.01 ft ³
Case Weight:	14.4 lbs
Cases per Skid:	40
Case GTIN#	10852482006554
Part GTIN#	00852482006557

Case measurements have a margin of error up to +/- .5 and are rounded to up to the nearest hundredth



Product dimensions listed are in inches and rounded to the nearest tenths

Manufactured by Oak Ridge Products, LLC
FDA Manufacturing Facility Registration
#1419823 & 3006060667

Oak Ridge Products, LLC.
ALL RIGHTS RESERVED

US:888-650-7444 Int'l:001-815-363-4700
www.oakridgeproducts.com

Rev. 10/26/2016