



5 Quart M-Style Sharps Container

Product Code: 0354-110M CLEAR

Product Code: 0354-150M Translucent RED

General Description



- Designed for in room use
- Suitable for unsupervised patient areas
- Meets all FDA, OSHA, and NOSHA regulations
- Fit 5 Qt cabinets and wall brackets
- Counter Balanced door allows for one hand use
- Large horizontal drop zone
- Locks for final disposal
- Use instructions on part label
- Made from High impact Polypropylene
- Comparable to: Bemis 175030, Covidien 8506SA, and 8507SA
- Fits Bemis 1051, and Covidien: 8515, 8518X, 8519C, 8516H, 8556H brackets

Packing Details

Part Dimensions:	10 x 4.4 x 9.2
Case:	20 each/case
Packaged:	Bulk, Nested
Case Dimensions:	25.50 x 17.50 x 11.00
Case Volume:	2.84 ft ³
Case Weight:	19 lbs
Cases per Skid:	30
Case GTIN#	10865273000082
Part GTIN#	00865273000085

Case measurements have a margin of error up to +/- .5 and are rounded to up to the nearest hundredth

For final lock: close door & press down on snaps



Product dimensions listed are in inches and rounded to the nearest tenths

Manufactured by Oak Ridge Products, LLC
FDA Manufacturing Facility Registration
#1419823 & 3006060667

Oak Ridge Products, LLC.
ALL RIGHTS RESERVED

US:888-650-7444 Int'l:001-815-363-4700
www.oakridgeproducts.com

Rev. 10/26/2016