



## STAT KIT® 750 Series

### EMERGENCY MEDICAL KITS

With essential and common ACLS medications and devices

Portable, pre-assembled emergency medical kits containing essential and common advanced cardiac life support medications<sup>1</sup> and devices in a well-organized case, so you are ready to respond to common medical emergencies.



# STAT KIT® 750 Series highlights

The STAT KIT 750 Series Emergency Medical Kits are designed to provide a comprehensive collection of medications and medical devices in a durable, mobile hard case. The contents are arranged into clearly labeled sections and medication trays, categorized by key areas of care for easy access during emergency situations.

- Essential and common advanced cardiac life support medications<sup>1</sup> and devices
- STAT KIT 750 version includes an adult epinephrine auto-injector
- STAT KIT 750-AI version includes adult and pediatric epinephrine auto-injectors
- Mobile, secure, watertight case with built-in wheels: 24.5" x 19.5" x 12.22", 40 pounds



## Always up to date

HealthFirst's Auto-Replenishment Service takes the guesswork out of staying up to date. The service uses smart automation to track the expiration dates of medications and devices and automatically replenishes them before expiration. Your practice benefits by saving time and gaining peace of mind.

## OnTraq™

OnTraq is an online account management tool that dynamically shows the status of your emergency medical kit and medications. It allows users to view each medication, expiration dates, lot numbers, order backorder status, and upcoming medication replenishment shipments. Rely on OnTraq's smart automation for dashboards, reports and task management.

## The STAT KIT 750 contains the following essential medications and devices:

### Medications

#### Top Tray:

Adenosine 3mg/mL, 2mL vial (5)  
Albuterol inhaler, 90mcg  
Amiodarone (50mg/ML) vial (3)  
Aspirin tablets 325mg (4)  
Atropine sulfate 0.1mg/mL prefilled syringe, 10mL (2)  
Dextrose 25% (Infant) 250mg/mL prefilled syringe, 10mL (1)  
Dextrose 50% 500mg/mL prefilled syringe, 50mL (1)  
Diphenhydramine 25mg capsule (2)  
Diphenhydramine 50mg/mL, 1mL vial (2)  
Dopamine 40mg/mL, 5mL vial (2)  
Epinephrine 1:1,000 1mg/mL, 1mL (2)  
Epinephrine 1:10,000 prefilled syringe, 10mL (5)  
Epinephrine auto-injector (Adult) 0.30mg (1)  
Epinephrine auto-injector (Pedi) 0.15mg (1)  
(STAT KIT 750-AI model)  
Flumazenil 0.1mg/mL, 10mL vial (1)  
Furosemide USP 100mg per 10mL, 10mL vial (2)  
Labetalol HCl USP 100mg/20mL, 20mL vial (1)  
Lidocaine HCl 1% and epinephrine 1:100,000, 20mL vial (1)  
Lidocaine 2% 20mg/mL prefilled syringe, 5mL (2)  
Magnesium sulfate USP 50% 5g/10mL prefilled syringe, 10mL (1)  
Naloxone 0.4mg/mL, 1mL vial (2)  
Nitroglycerin 0.4mg tablets (25)  
Ondansetron 2mg/mL, 2mL vial (2)  
Oral glucose gel, 15g  
Solu-Medrol® 125mg, 2mL Act-O-Vial® (1)  
Verapamil 2.5mg/mL, 2mL vial (2)

#### Bottom Tray:

Lactated Ringer's and 5% Dextrose 1000mL bag (1)  
Sodium chloride 0.9% 500mL bag (2)

### Syringes, Needles,

#### IV Equipment

Syringes with needle:  
3cc 21G x 1 1/2" (1)  
3cc 20G x 1 1/2" (1)  
10cc 21G x 1 1/2" (1)  
10cc 20G x 1 1/2" (1)  
Syringes without needle:  
1cc syringe (2)  
3mL syringe (1)  
12mL syringe (1)  
30mL syringe (1)  
Needles:  
25G x 5/8" (1)  
25G x 1 1/2" (1)  
22G x 1 1/2" (1)  
20G x 1 1/2" (2)  
18G x 1 1/2" (1)  
18G x 1 1/2" Blunt filter (2)  
IV catheter 22G x 1" (1)  
IV catheter 20G x 1" (2)  
IV administration tubing  
Alcohol wipes (6)  
Tape, hypoallergenic (1 roll)  
Tourniquet, latex free (2)  
  
1 quart sharps disposal container

### Airway Equipment

Laryngeal Mask Airways (LMA)  
Size 3 (1)  
Size 4 (1)  
Size 5 (1)  
Nasopharyngeal airway set (one each):  
26FR, 28FR, 30FR, 32FR, 34FR  
Oropharyngeal airways:  
60mm (1)  
80mm (1)  
100mm (1)  
Laryngoscope (1)  
Size 2 Macintosh blade (1)  
Size 4 Macintosh blade (1)  
Battery, C cell (1 pair)  
Endotracheal tube set  
3mm, uncuffed (1)  
5mm, cuffed (1)  
7mm, cuffed (1)  
Hand-held suction unit (1)  
Yankauer suction tip (1)  
Magill forceps (1)

### Resuscitation & Oxygen

Bag-Valve-Mask (BVM) adult resuscitator (1)  
Adult mask (1)  
Toddler mask (1)  
CO<sub>2</sub> detector (Adult)

### Wound Care

Gauze pads, 4 x 4" (4)  
Gauze sponges, 3 x 3" (4)  
Gauze bandage, 3" roll (2)  
Compression bandage, 3" roll (2)  
Chloraprep® swabstick, 1.75mL (1)  
Betadine prep pads (2)

### Instruments

Scalpel (1)  
Shears (1)

### Vitals, PPE, Guides

Stethoscope (1)  
Sphygmomanometer (1)  
Thermometer strips (2)  
Headlamp (1)  
Nitrile gloves, pair (6)  
N95 respirator (1)  
CPR shield (1)  
Safety glasses (1)  
ACLS pocket guide

### Miscellaneous

Notepad & pen  
Security seals (4)

<sup>1</sup> American Heart Association 2025 Handbook of Emergency Cardiovascular Care for Healthcare Providers

A physician license is required for the purchase of prescription medications and devices.

IMPORTANT: Medications will be shipped in the available form factor at the time of fulfillment. Some medications or devices may not be available at the time of order due to national supply issues. Kits may ship with backordered items. Backordered items will ship as soon as they are available. Visit [healthfirst.com](http://healthfirst.com) for complete terms and conditions of sale.