

EDAN Monitor Comparison Chart

	Vital Signs Monitors		Patient Monitors with ECG and Capnography				
	M3 Item #1023650	iM3 Item #1023680	CX12 Item #1025010	iM60GT Item #1023050	iM70GT Item #1023070	X10 Item #1023060	X12 Item #1023430
<i>HealthFirst is an authorized distributor of the following EDAN monitors. This chart is for general comparison. Please refer to the manufacturer spec sheets for more details.</i>							
Blood Pressure (NIBP)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Pulse Rate (PR)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Pulse Oximetry (SpO ₂)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
ECG 3/5 Lead	No	No	Yes	Yes	Yes	Yes	Yes
Arrhythmia and ST segment analysis	No	No	Yes	Yes	Yes	Yes	Yes
Built-in Capnography (EtCO ₂)	No	No	Yes	Yes	Yes	Yes	Yes
Temperature	No	No	Yes	Yes	Yes	Yes	Yes
Screen Size	5.6" LCD	8.0" Touchscreen	13.3" with Quick Gesture Interface	10.4" Touchscreen	12.1" Touchscreen	10.1" Touchscreen	12.1" Touchscreen
Size/Dimensions	6.85"W x 9.25"H x 7.44"D	6.26"W x 10.31"H x 6.54"D	13"W x 9.61"H x 6.94"D	12.00"W x 10.00"H x 6.34"D	13.00"W x 11.22"H x 6.22"D	10.28"W x 9.69"H x 5.75"D	13.15"W x 10.40"H x 6.38"D
Weight	7.72 pounds	5.51 pounds	8.82 lbs	8.38 pounds	10.36 pounds	6.17 pounds	9.14 pounds
User interface	Programmable keys	Customizable shortcuts	Customizable shortcuts	Customizable shortcuts	Customizable shortcuts	Customizable shortcuts	Customizable shortcuts
Printer	Integrated thermal printer	Integrated thermal printer	Integrated thermal printer	Integrated thermal printer	Integrated thermal printer	Integrated thermal printer	Integrated thermal printer
Printer Location	Built-in	Built-on	Built-in	Built-in	Built-in	Built-on	Built-on
Wi-Fi	Built in Wi-Fi	Built in Wi-Fi	Built in Wi-Fi	Built in Wi-Fi	Built in Wi-Fi	Built in Wi-Fi	Built in Wi-Fi
USB Ports	1 USB port	1 USB, 1 Micro USB port	2 USB Ports	2 USB ports	2 USB ports	2 USB ports	2 USB ports
Video Output	No	No	HDMI Video Output	VGA Output	VGA Output	VGA Output (Optional)	VGA Output (Optional)
Memory	Internal and external memory	Internal and external memory	Internal and external memory	Internal and external memory	Internal and external memory	Internal and external memory	Internal and external memory
NIBP Storage	Up to 30,000 text-based sets	1,200 sets	1,600 sets	1,200 sets	1,200 sets	1,200 sets	1,200 sets
Battery	Lithium-ion rechargeable battery	Lithium-ion rechargeable battery	Lithium-ion rechargeable battery	Lithium-ion rechargeable battery	Lithium-ion rechargeable battery	Lithium-ion rechargeable battery	Lithium-ion rechargeable battery
Connectors	Lemo	7-Pin	7-Pin	Lemo	Lemo	7-Pin	7-Pin
Warranty	2 years	2 years	3 years	3 years	3 years	3 years	3 years
Comes standard with: <i>*replacement items sold individually</i>	3 adult BP cuffs (large adult, adult, small adult), BP hose, SpO ₂ sensor, and printer.		3 adult BP cuffs (large adult, adult, small adult), BP hose, SpO ₂ sensor, ECG cable, electrodes, temperature probe, grounding cable, straight sampling line, sampling cannula, dual sampling cannula, water trap, printer and printer paper.				

The information contained in this sales sheet is provided for informational purposes only and is subject to change without notice. All products are sold as-is, and the distributor makes no warranties or representations, either express or implied, regarding the accuracy, completeness, or suitability of the products for any particular purpose. All prices and product availability are subject to change. It is the responsibility of the purchaser to verify product specifications and ensure the products are suitable for their intended use.

